

IHC PROTOCOL - IMMUNOFLUORESCENCE DETECTION (IHC-IF)

Protocol for immunohistochemistry with fluorescent detection optimized for Triple A Polyclonals™ and PreciSA Monoclonals™ from Atlas Antibodies.

MATERIALS NEEDED

Cryosections mounted on SuperFrost slides

- Target retrieval solutions pH6 or pH9 (e.g. DAKO or ThermoFisher)
- Incubation chamber with wet Wettex stripes
- DAKO Cytomation pen
- 1xPBS
- Fluorophore-conjugated secondary antibodies for regular immunofluorescence
- Mounting media (e.g. ProLong Gold with or without DAPI)
- Coverslips
- Eppendorf tubes for dilution of antibodies
- Pipettes and pipette-tips.

CRYO SECTIONS

- **Perfusion-fixed**, 10-30% sucrose cryoprotected sections cut in a microtome at 14 µm, thaw-mounted on Super Frost slides. Dry sections on slides for additional 2 h at room temperature (RT). Rehydrate in 1xPBS 2x15 min, or proceed with antigen retrieval.
- **Snap-frozen** sections cut in a microtome at 14 µm thaw-mounted on Super-Frost slides and dried for 2h. Pre-fix sections in ice-cold paraformaldehyde (4% PFA) for 20 min (put frozen sections directly into 4% PFA, do not allow thawing). Rinse in 1x PBS 15 min.
-

ANTIGEN RETRIEVAL FOR CRYO SECTIONS (OPTIONAL)

If target requires antigen retrieval, standard target retrieval (HIER) solutions can be used, but the temperature and time should be reduced.

HEAT-INDUCED EPI TOPE RETRIEVAL

1. DAKO or Thermo HIER solution pH 6, a modified citrate buffer, 60-70°C, 10 min.
2. DAKO or Thermo HIER solution pH 9, a Tris/EDTA buffer, 60-70°C, 10 min.

NOTE: The specified working dilutions of the primary antibodies are to be considered as a guideline only. Optimal dilutions must be determined by the user.

IHC PROCEDURE

DAY I

- Draw a circle around each section with DAKO Cytomation pen.
- Rinse slides in 1x PBS for 15 min
- Block sections in 2% normal serum of the secondary antibody host, in PBS 30 min at RT (optional)
- Add primary antibodies (diluted in 1x PBS, containing 0.3% Triton, 0.01% NaAzide, 0.02% Bacitracin), approx. 180µl/slide.
- Incubate in a humidified chamber overnight at 4°C.

DAY II

- Rinse slides in 1x PBS for 15 min at RT (under a black cover)
- Incubate with secondary antibody conjugated with fluorophore (1:80 -1:400 in PBS), 30 min at 37°C or 1-2 h at RT.
- Wash in PBS 2x15 min at RT under a black cover
- Mount in anti-fading solution (e.g ProLong Gold)
- Leave slides for 24h at 4C before microscopy/scanning
- Store at 4C

Warranty: The products supplied by Atlas Antibodies are warranted to meet stated product specifications and to confirm label descriptions when used and stored properly. Unless otherwise stated, this warranty is limited to one year from the date of the sale for products used, handled, and stored according to Atlas Antibodies AB's instructions. Atlas Antibodies AB's sole liability is limited to the product replacement or refund of the purchase price. All products are supplied for research use only. They are not intended for medicinal, diagnostic, or therapeutic use. No products from Atlas Antibodies may be resold, modified for resale, or used to manufacture commercial products without prior written approval from Atlas Antibodies AB.

Atlas Antibodies logo, Triple A Polyclonals, PreciSA Monoclonals, and PrEST Antigens are trademarks or registered trademarks of Atlas Antibodies AB. All other trademarks are the property of their respective owners. Products are for research use only. Not for use in diagnostic procedures. © Atlas Antibodies AB 2022.