

PRODUCT SPECIFICATION

Product Name	Anti-NHS
Product Number	HPA031497
Gene Description	Nance-Horan syndrome (congenital cataracts and dental anomalies)
Clonality	Polyclonal
Isotype	IgG
Host	Rabbit
Antigen Sequence	TDLPYLEESTLTTAALSPSKIRPHTANKSVSRQYSTEDTILSFLDSSAVE MGPDKLHLEKNSTFDVKNRCDPETITSAGSSLLDSNVTKDQVRTETETPIP ENTPTKNCAFPTEGFQRVSAARPNDLDGKIIQYGPDPDETLEQ
Purification Method	Affinity purified using the PrEST antigen as affinity ligand
Verified Species Reactivity	Human
Recommended Applications	IHC (Immunohistochemistry) - Antibody dilution: 1:20 - 1:50 - Retrieval method: HIER pH6
Characterization Data	Available at atlasantibodies.com/products/HPA031497
Buffer	40% glycerol and PBS (pH 7.2). 0.02% sodium azide is added as preservative.
Concentration	Lot dependent
Storage	Store at +4°C for short term storage. Long time storage is recommended at -20°C.
Notes	Gently mix before use. Optimal concentrations and conditions for each application should be determined by the user.

 *Product of Sweden. For research use only. Not intended for pharmaceutical development, diagnostic, therapeutic or any in vivo use. No products from Atlas Antibodies may be resold, modified for resale or used to manufacture commercial products without prior written approval from Atlas Antibodies AB.*

Warranty: The products supplied by Atlas Antibodies are warranted to meet stated product specifications and to conform to label descriptions when used and stored properly. Unless otherwise stated, this warranty is limited to one year from date of sales for products used, handled and stored according to Atlas Antibodies AB's instructions. Atlas Antibodies AB's sole liability is limited to replacement of the product or refund of the purchase price. All products are supplied for research use only. They are not intended for medicinal, diagnostic or therapeutic use. No products from Atlas Antibodies may be resold, modified for resale or used to manufacture commercial products without prior written approval from Atlas Antibodies AB Rev. December 2012.